

Leadgene-XenoForm Tumor Gel

- XenoForm Tumor Gel offers a secure and stable extracellular matrix microenvironment. It is constructed from high-activity synthetic biological materials and has excellent biocompatibility and an adjustable density. The gel, which is convenient to handle at room temperature, stands out as the optimal choice for conducting tumor growth experiments in mouse models.

Characteristics



► Bio-Synthesis



► Stable



► Tunable



► R.T.




► Easy

Comparison

Products	XenoForm Tumor Gel	Comparison-1	Comparison-2
Origin	Synthetic biological material	Basement membrane extracted from the EHS mouse sarcoma	Animal-free hydrogel
Growth Factor	Without	Complicated	Without
Gel density	Adjustable	-	Adjustable
Gelation	Room temperature	>10°C	Room temperature
Gelation time	10 minutes	-	20-30 minutes
Temperature before injection	Room temperature	4°C	Room temperature
Tumor formation	4 weeks	4 weeks	6-8 weeks

Procedure



Thoroughly mix the gel with the Cell/R buffer in a 1:1 ratio.




Draw the gel mixture into a syringe and allow it to react at room temperature for at least 10 minutes.



Inject the required volume into the subcutaneous tissue the experimental mice (0.2 mL).



Harvest tumor tissue after 4-5 weeks.



Apply to in vivo tumorigenesis tests in various cancer models.