

## M-MLV Reverse Transcriptase (Glycerol-Free)

Cat no. LDG0011RF

### Product Overview

#### Package component

Specification	Item	Amount
20,000 U	M-MLV Reverse Transcriptase (Glycerol-Free)	1 vial (200 U/μL)
50,000 U	M-MLV Reverse Transcriptase (Glycerol-Free)	1 vial (200 U/μL)

### Description

Moloney Murine Leukemia Virus (M-MLV) Reverse Transcriptase is an RNA-dependent DNA polymerase that synthesizes the first strand of complementary cDNA from a single-stranded RNA template with hybridized primer. This kit features high activity formulation of M-MLV RT and enzyme dilution buffer. The enzyme formulation does not contain glycerol and is compatible for further lyophilization process.

### Source

*Escherichia coli*

### Activity

One unit of M-MLV Reverse Transcriptase is defined as the amount of the enzyme incorporates 1 nmol of dTTP into acid-insoluble product in 10 minutes at 37°C.

### Storage and Stability

Stored at -20°C. Avoid repeated freeze/thaw cycles.

### Procedure

#### First strand cDNA synthesis:

The following procedure is a general guideline for RT-qPCR reaction. To maintain an RNase-free environment, always wear disposable gloves, and use laboratory consumables

and water of nuclease-free grade during the whole experiment course.

#### RT-qPCR reaction set-up:

- Place all required reagents on ice.

Component	Amount	Final concentration
2X Probe qPCR Master Mix	10 μL	1X
M-MLV Reverse Transcriptase (Glycerol-Free)	1 μL	200 U/rxn
Forward primer (10 μM)	0.8 μL	0.4 μM
Reverse primer (10 μM)	0.8 μL	0.4 μM
Probe (10 μM)	0.4 μL	0.2 μM
RNA template	X μL	≤ 1 μg (total RNA)
Nuclease-Free H <sub>2</sub> O	Y μL	-
<b>Total reaction volume</b>	20 μL	-

- Gently mix the reaction thoroughly to achieve uniform distribution and briefly centrifuge.
- Thermal cycling conditions for standard qPCR

Step	Cycles	Temperature	Time
Reverse transcription	1	42 – 50°C	10 – 15 min
Enzyme activation	1	95°C	5 min
Denaturation	40 - 45	95°C	5 – 15 sec
Annealing/Extension		55 – 65 °C	30 – 60 sec

### Important notes

After the reaction is complete, M-MLV RTase can be inactivated by incubation at 65°C for 20 minutes.

#### Primer/Probe concentration

Final concentrations of 400 nM (each primer) are suitable for most reactions. To obtain optimal condition, primer

concentration can be titrated between 0.2 – 1  $\mu$ M.

A final concentration of 200 nM (probe) is suitable for most reactions. To obtain optimal condition, probe concentration can be titrated between 0.1 – 0.3  $\mu$ M.

#### **Annealing/Extension optimization**

To obtain optimal condition, annealing/extension temperature can be adjusted between 55°C – 65°C, annealing/extension time can be extended up to 60 sec.

#### **Target length**

Appropriate amplicon length should be arranged between 80 – 200 bp.

#### **Disclaimer**

This product is for research use only and is not intended for diagnostic use.

## **LEADGENE BIOMEDICAL, INC.**

No.9, Ln. 147, Zhengbei 1st Rd., Yongkang Dist., Tainan City 710, Taiwan R.O.C. TEL: +886-6-2536677 FAX: +886-6-2531536

[www.leadgenebio.com](http://www.leadgenebio.com)