

## Bst DNA Polymerase (Large Fragment) (Glycerol-Free)

Cat no. LDG0012RF

### Product Overview

#### Package component

Specification	Item	Amount
1,600 U	Bst DNA Polymerase (Large Fragment) (Glycerol-Free)	1 vial (8 U/μL)
	10X Bst DNA Polymerase Reaction Buffer (Mg <sup>2+</sup> free)	1 vial (1 mL)
	100 mM MgSO <sub>4</sub>	1 vial (1 mL)
8,000 U	Bst DNA Polymerase (Large Fragment) (Glycerol-Free)	1 vial (8 U/μL)
	10X Bst DNA Polymerase Reaction Buffer (Mg <sup>2+</sup> free)	3 vial (1 mL)
	100 mM MgSO <sub>4</sub>	3 vial (1 mL)

### Description

Bst DNA Polymerase (Large fragment) is an enzyme of *Bacillus stearothermophilus* DNA polymerase which can catalyze 5' → 3' polymerase activity but lacks 5' → 3' exonuclease activity. Bst DNA Polymerase offers strand displacement capabilities, making it ideal for isothermal amplification.

The enzyme formulation does not contain glycerol and is compatible for further lyophilization process.

### Source

*Escherichia coli*

### Activity

One unit of Bst DNA Polymerase is defined as the amount of the enzyme incorporates 10 nmol of dNTP into acid-insoluble product in 30 minutes at 65°C.

### Storage and Stability

Stored at -20°C. Avoid repeated freeze/thaw cycles.

### Procedure

LAMP reaction recipe:

- Place all required reagents **on ice** and add each of them following the order suggested below.

Component	Amount	Final concentration
10X Bst DNA Polymerase Reaction Buffer (Mg <sup>2+</sup> free)	2.5 μL	1X
100 mM MgSO <sub>4</sub>	0.5-2.5 μL	2-10 mM
10 mM dNTP mix	3.5 μL	1.4 mM each
10X FIP/BIP primers	1 μL	1.6 μM
10X F3/B3 primers	1 μL	0.2 μM
10X LoopF/B primers	1 μL	0.8 μM
DNA template	X μL	10 copies or more
Nuclease-Free H <sub>2</sub> O	Y μL	-
Bst DNA Polymerase (Large Fraction) (Glycerol-Free) (8 U/μL)	1 μL	8 U/rxn
<b>Total reaction volume</b>	<b>25 μL</b>	<b>-</b>

- Gently mix the reaction thoroughly to achieve uniform distribution.
- Incubate at 65°C for 30-60 minutes.
- MgSO<sub>4</sub> (2-10 mM), Bst DNA Polymerase (40-320 U/mL) and temperature (50-65 °C) can be adjusted for optimal results.
- Reaction preparations may be scaled up or down

proportionately.

### Important notes

It is not recommended to perform reaction above 70 °C. Bst DNA Polymerase cannot be used for thermal cycle sequencing.

### Disclaimer

This product is for research use only and is not intended for diagnostic use.

## LEADGENE BIOMEDICAL, INC.

No.9, Ln. 147, Zhengbei 1st Rd., Yongkang Dist., Tainan City 710, Taiwan R.O.C. TEL: +886-6-2536677 FAX: +886-6-2531536  
[www.leadgenebio.com](http://www.leadgenebio.com)