



Product Information & Manual

Information of other products is available at: www.leadgenebio.com

T7 RNA Polymerase (200 U/μL)

Cat no. LDG0001RI

Product Overview

Package component

Item	Amount
T7 RNA Polymerase (200 U/μL)	25,000 U
10X RNA Polymerase reaction buffer	1 vial (1 mL)
100 mM DTT	1 vial (1 mL)

Description

Bacteriophage T7 RNA Polymerase is a DNA-dependent RNA polymerase with high specificity for the T7 promoter. This enzyme catalyzes the 5' → 3' synthesis of RNA from DNA downstream from its promoter.

Source

Escherichia coli

Activity

One unit of T7 RNA Polymerase is defined as the amount of the enzyme incorporates 1 nmol of ATP into acid-insoluble product in 1 hour at 37°C.

Storage buffer

The enzyme is supplied in 100 mM Tris-HCl (pH 7.9), 20 mM KCl, 1 mM DTT, 1 mM EDTA, 0.1% Triton® X-100 and 50% (v/v) glycerol.

10X RNA Polymerase Reaction Buffer

400 mM Tris-HCl (pH 8.0), 60 mM MgCl₂, and 20 mM spermidine.

Storage and Stability

Stored at -20°C. For optimal storage, aliquot the enzyme, reaction buffer and DTT reagent into smaller quantities and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple

freeze/thaw cycles.

Procedure

Standard RNA synthesis procedures:

- Below reaction mixture should be prepared under room temperature and combined in the following order:

Component	Amount	Final concentration
Nuclease-Free H ₂ O	X μL	-
Template DNA	0.5-1 μg	-
10X RNA Polymerase Reaction Buffer	2 μL	1X
ATP (100 mM)	0.6 μL	3 mM
UTP (100 mM)	0.6 μL	3 mM
CTP (100 mM)	0.6 μL	3 mM
GTP (100 mM)	0.6 μL	3 mM
100 mM DTT	2 μL	10 mM
T7 RNA Polymerase (200 U/μL)	1 μL	-
RNase inhibitor (optional)	0.5 μL	1 U/μL
Total reaction volume	20 μL	-

- Incubate at 37°C for 30 minutes to 2 hours.
- Above reaction mixture may be scaled up or down proportionately.

Important notes

- Transcription reaction should be performed under RNase free condition. Use nuclease-free tubes, reagents, and water to avoid RNase contamination. Also, wear gloves when working with RNA.
- To obtain optimal condition, NTP concentration can be titrated between 3 – 5 mM.
- The volume of T7 RNA Polymerase can be titrated between 1-2 μL in the IVT reaction to optimize your assay.

Disclaimer

This product is for research use only and is not intended for diagnostic use.

LEADGENE BIOMEDICAL, INC.

No.9, Ln. 147, Zhengbei 1st Rd., Yongkang Dist., Tainan City 710, Taiwan R.O.C. TEL: +886-6-2536677 FAX: +886-6-2531536
www.leadgenebio.com