



## Product Information & Manual

Information of other products is available at: [www.leadgenebio.com](http://www.leadgenebio.com)

### LEADSPHERE® One-Step RT-PCR Master Mix

Cat no. LDG0027RF

#### Product Overview

##### Package component

Item	Content
LEADSPHERE® One-Step RT-PCR Master Mix	Inquiry

#### Description

Essential components for reliable RT-PCR are included in each single LEADSPHERE® One-Step RT-PCR Master Mix, which saves time for reaction set-up and reduces the risk of contamination and pipetting errors. In addition, the stabilized format facilitates shipment and storage without cooling. LEADSPHERE® One-Step RT-PCR Master Mix includes reverse transcriptase, DNA polymerase, high-quality dNTPs, and optimized PCR buffer. This product is suitable for probe-based detection and endpoint DNA amplification. The reaction can be set up by reconstituting the LEADSPHERE™ One-Step RT-PCR Master Mix by adding PCR-grade water along with the RNA template, primers, and probe to a total volume of 20 µL.

#### Storage and Stability

Stored at 4-30°C and keep dry.

#### Procedure

The following procedure is a general guideline for a One-step RT-qPCR reaction. To maintain an RNase-free environment, always wear disposable gloves, and use laboratory consumables and water of nuclease-free grade during the whole experiment course.

RT-qPCR reaction set-up:

- Place all required reagents on ice

Component	Amount	Final concentration
LEADSPHERE® One-Step RT-PCR Master Mix	1 tablet	1 X
Forward primer (10 µM)	0.8 µL	0.4 µM
Reverse primer (10 µM)	0.8 µL	0.4 µM
Probe (10 µM)	0.4 µL	0.2 µM
RNA template	X µL	≤ 1 µg (total RNA)
Nuclease-Free H <sub>2</sub> O	Y µL	-
<b>Total reaction volume</b>	<b>20 µL</b>	<b>-</b>

- Prepare the above reagents and gently mix to 20 µL total reaction volume to achieve uniform distribution and briefly centrifuge.
- Thermal cycling conditions for standard qPCR

Step	Cycles	Temperature	Time
Reverse transcription	1	50°C	10 – 15 min
Enzyme activation	1	95°C	5 min
Denaturation	40-50	95°C	5 – 15 sec
Annealing/Extension		55-65°C	30 – 60 sec

#### Important notes

##### Primer/Probe concentration

Final concentrations of 400 nM (each primer) are suitable for most reactions. To obtain optimal condition, primer concentration can be titrated between 0.2-1 µM.

A final concentration of 200 nM (probe) is suitable for most reactions. To obtain optimal condition, probe concentration can be titrated between 0.1-0.3 µM.

##### Annealing/Extension optimization

To obtain optimal condition, annealing/extension temperature can be adjusted between 55°C-65°C, annealing/extension time can be extended up to 60 sec.

##### Target length

Appropriate amplicon length should be arranged between 80-200 bp

## Disclaimer

This product is for research use only and is not intended for diagnostic use.

## LEADGENE BIOMEDICAL, INC.

No.9, Ln. 147, Zhengbei 1st Rd., Yongkang Dist., Tainan City 710, Taiwan R.O.C. TEL: +886-6-2536677 FAX: +886-6-2531536  
[www.leadgenebio.com](http://www.leadgenebio.com)