



## Product Information & Manual

Information of other products is available at: [www.leadgenebio.com](http://www.leadgenebio.com)

### HyLink™ Biotin Labeling Kit (SpinDesalt Column), 100 µg\*10

Cat no. LDG0007RC

#### Product Overview

| Components                       |                |
|----------------------------------|----------------|
| Biotin                           | 1 mg, 1 vial   |
| Biotin reconstitution buffer     | 300 µL, 1 vial |
| 10X Modifier                     | 300 µL, 1 vial |
| PBS                              | 10 mL, 1 vial  |
| SpinDesalt Column<br>(LDG0008RC) | 10 vials       |

#### Description

Biotin is a widely used and powerful tool for research due to the specific and high affinity with streptavidin/avidin. Antibody or protein conjugated with several biotin molecules can also amplify the detection signal through streptavidin-conjugated molecule such as streptavidin-HRP, streptavidin-FITC, etc. Biotinylated antibody or protein can be used in various applications including ELISA, WB, IHC, IFA and FACS. Leadgene HyLink™ Biotin Labeling Kit (SpinDesalt Column) is designed for biotinylation of antibody or protein (100 µg-1 mg). It provides a rapid and easy process with high efficiency to conjugate biotin to antibody or protein.

#### Materials required but not provided

- (1) Disposable microcentrifuge tubes (1.5 or 2 mL).
- (2) 5 µL to 1000 µL adjustable single-channel micropipettes with disposable tips.
- (3) Timer
- (4) Vortex mixer

- (5) Incubator capable of maintaining temperature at 37±1°C.

#### Storage and Stability

- (1) Stored at -20°C except for SpinDesalt Column stored at 4°C. Avoid repeated freeze/thaw cycles.
- (2) Equilibrate kit to room temperature before use.
- (3) The kit is stable for one year under proper storage conditions.

#### Procedure

##### Biotin conjugation protocol

##### ● For antibody conjugation

- (1) Dissolve the 1 mg Biotin with 200 µL Biotin reconstitution buffer to a concentration of 5 mg/mL.
- (2) Dissolve antibody in PBS (**Refer to Important notes point 1 before dissolving**) or other buffers that do not contain amine, Tris or glycerol. Use 10X Modifier (e.g. Add 1 µL of 10X Modifier for 9 µL of antibody) to adjust the antibody buffer condition.
- (3) Make sure all buffers are well dissolved before using. If not, please vortex the vial to make salts dissolved.
- (4) Add Biotin into the mixture that containing antibody and 10X Modifier (e.g. **Add 20 µL of Biotin for 100 µg of antibody**). Mix gently by pipetting several times.
- (5) Next, incubate in the dark at room temperature for 1 hour, then move to 4°C for 2 hours.
- (6) After incubating, using SpinDesalt Column to remove excess biotin (**Refer to Important notes point 2 before using column**). The conjugates can be used after desalting.

- For protein conjugation, the amount of protein can be calculated by formula below:

Quantities of protein =

$$\text{Biotin conjugation (e.g. 1 mg)} \times \frac{\text{M.W. (KDa) protein}}{150 \text{ kDa IgG}}$$

## LEADGENE BIOMEDICAL, INC.

No.9, Ln. 147, Zhengbei 1st Rd., Yongkang Dist., Tainan City 710, Taiwan R.O.C. TEL: +886-6-2536677 FAX: +886-6-2531536  
www.leadgenebio.com

### SpinDesalt Column protocol

- (1) Prepare a **SpinDesalt Column** by breaking off the bottom closure and placing the column into a microcentrifuge tube.
- (2) Centrifuge the column at  $1,000 \times g$  for 1 minute, discard the storage buffer and return column to the same microcentrifuge tubes.
- (3) Adding 0.25 mL of **PBS** to the top of the resin bed and centrifuging at  $1,000 \times g$  for 1 minute. Discard the flowthrough and repeat this step 3 times.
- (4) Place the column into a new microcentrifuge tube and apply approximately 0.1-0.25 mL of the conjugates directly onto the resin bed. Centrifuge the column at  $1,000 \times g$  for 1 minute.
- (5) The collected flowthrough solution is the purified conjugates.

### 2. SpinDesalt Column

- (1) For conjugates < 0.1 mL, please use PBS to adjust the volume to at least 0.1 mL to increase the recovery rate of SpinDesalt Column.
- (2) The recovery rate of the SpinDesalt Column is related to the type of protein and other biomolecules, usually exceeding 85%. Increasing the sample concentration or volume can improve the recovery rate.
- (3) After completing the conjugation, the resin bed in the column can be temporarily stored in PBS.

### Disclaimer

This product is for research use only and is not intended for diagnostic use.

### Important notes

#### 1. Biotin labeling

- (1) Antibody concentrations of 1-2 mg/mL generally give optimal results. Recommended amount and volume of antibody for optimal results.

| Kit size    | Antibody amount | Reaction volume |
|-------------|-----------------|-----------------|
| 100 µg x 10 | 0.1 mg-2 mg     | 100-250 µL      |

- (2) Common non-buffering salts (e.g. sodium chloride) have no effect on conjugation efficiency. Avoid buffer component that contains primary amine (e.g. amino acid or ethanolamine) and thiols (e.g. 2-Mercaptoethanol or DTT).

Components that have on effect or little effect on labeling reaction:

- up to 50 mM Tris
- up to 50 mM HEPES
- up to 10% glycerol
- up to 0.02% sodium azide