

## HyLink™ Biotin Labeling Kit, 100 µg\*3

**Catalog Number** LDG0005RC

**Package** 100 µg\*3 / Customized package

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

Biotin is a widely used and powerful tool for research due to the specific and high affinity with streptavidin. Antibody or protein conjugated with several biotin molecules can also amplify the detection signal through streptavidin-conjugated molecule such as streptavidin-HRP, streptavidin-FITC, etc. Biotinylated antibody or protein can be used in various applications including ELISA, WB, IHC, IFA and FACS. Leadgene HyLink™ Biotin Labeling Kit is designed for biotinylation of a small quantity (10 µg-1 mg) of antibody or protein. It provides a rapid and easy process with high efficiency to conjugate biotin to antibody or protein.

#### Components

| Package    | Items        | Quantity |
|------------|--------------|----------|
| 100 µg × 3 | Biotin       | 3 vials  |
|            | 10× Modifier | 1 vial   |
|            | 10× Quencher | 1 vial   |

### Instruction

#### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

#### Stability & Storage

This product is stable after storage at:

- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Equilibrate kit to room temperature before use.

#### Tainan Headquarter

☎ +886-6-2536677 📠 +886-6-2531536

✉ [info@leadgene.com.tw](mailto:info@leadgene.com.tw)

#### Taipei Office

☎ +886-2-27065528

**Disclaimer :** For Research Use or Further Manufacturing Only.

**Tainan Headquarter**

+886-6-2536677 +886-6-2531536

info@leadgene.com.tw

**Taipei Office**

+886-2-27065528