

NTP Set, 100 mM Solutions

Catalog Number LDG0006RI

Package 1 mL*4 / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

The NTP Set consists of 4 separate 100 mM solutions of ATP, UTP, CTP, and GTP (pH 7.5) as sodium salts.

Product Note

Mix thoroughly before use and dilute with nuclease-free water if necessary.

Specifications

Application

In vitro transcription, RNA amplification, miRNA and siRNA synthesis

Concentration

100 mM

Purity

≥98% purity of each component confirmed by HPLC.

Form

Liquid. Clear aqueous solutions (pH 7.4).

Instruction

Tainan Headquarter

+886-6-2536677 +886-6-2531536

info@leadgene.com.tw

Taipei Office

+886-2-27065528

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

Stability & Storage

This product is stable after storage at:

- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarter

☎ +886-6-2536677 📠 +886-6-2531536

✉ info@leadgene.com.tw

Taipei Office

☎ +886-2-27065528