

## Anti-EGFP Antibody Sepharose Purification Kit

<b>Catalog Number</b>	LDG0018RD
<b>Package</b>	5 rxn / 10 rxn / Customized package

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

Anti-EGFP sepharose, is consist of an EGFP VHH coupled to NHS-sepharose. It is an efficient technique for isolating recombinant proteins or mammalian expression proteins. The EGFP epitope system relies on a VHH EGFP recombinant antibody, which is able to react with N- and C- terminal GFP tagged and EGFP tagged fusion proteins. This product can be used for the immunoprecipitation or immune affinity purification.

The purified antibody is immobilized at 3-5 mg antibody per mL of 50% slurry and this kit allows a rapid and efficient affinity purification of active EGFP fusion proteins. The affinity resin allows an efficient binding of EGFP fusion proteins without the need for preliminary steps and calibrations. The affinity bound EGFP fusion proteins can be efficiently eluted from the resin by acid condition. The eluted proteins can be used for characteristic analysis.

#### Product Note

The optimal protein concentration of the cell lysate should be determined by the end user.

### Specifications

#### Clonality

VHH

#### Application

Immunoprecipitation, Immunoaffinity Purification

#### Conjugation

Sepharose

#### Specificity

Recognizes GFP tag and EGFP tag on either the N- or C-terminus of the recombinant fusion protein.

#### Sample Type

Cell lysate

#### Binding capacity

The binding and elution capacity of 1 mL settled Anti-EGFP sepharose may be vary, commonly more than 1 mg GFP or EGFP fusion proteins.

#### Tainan Headquarters

+886-6-2536677

bd@leadgene.com.tw

#### Innovation & Research Center

+886-2-27065528

#### CLD Center

+886-6-2536677

## Instruction

### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C.

### Stability & Storage

This product is stable after storage at:

- 2-8°C for unopened product

Please refer to product manual for storage constructions

**Disclaimer :** For Research Use or Further Manufacturing Only.

### Tainan Headquarters

+886-6-2536677

✉ [bd@leadgene.com.tw](mailto:bd@leadgene.com.tw)

### Innovation & Research Center

+886-2-27065528

### CLD Center

+886-6-2536677