

# **Annexin V-FAM**

Catalog Number	LDG0004RB
Package	25 rxn / 50 rxn / 100 rxn / Customized package

For full product information, images and publications, please visit our website.



## **Overview**

## **Description**

Leadgene fluorescent dye (FAM) conjugated Annexin V is highly purified product. During early apoptosis, cells will translocate membrane phosphatidylserine (PS) from the inner face of the membrane to the cell surface. The product can be used in one-step staining procedure without wash step within 20 minutes.

#### **Components**

Package	Items	Quantity
25 rxn	Annexin V-FAM	1 vial (0.125 mL)
50 rxn	Annexin V-FAM	1 vial (0.25 mL)
100 rxn	Annexin V-FAM	1 vial (0.5 mL)

# **Specifications**

**Application** Reactivity Mammalian Cell apoptosis assay **Detection Method Sample Type** Fluorescent Adherent cells | Suspension cells

## Instruction

Tainan Headquarter

**Taipei Office** 



## **Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

## Stability & Storage

This product is stable after storage at:

• 2-8°C for 12 months under sterile conditions from date of receipt.

Protect from light.

**Disclaimer:** For Research Use or Further Manufacturing Only.