

Annexin V-FAM

Catalog Number	LDG0004RB
Package	25 rxn / 50 rxn / 100 rxn / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

Leadgene fluorescent dye (FAM) conjugated Annexin V is highly purified product. During early apoptosis, cells will translocate membrane phosphatidylserine (PS) from the inner face of the membrane to the cell surface. The product can be used in one-step staining procedure without wash step within 20 minutes.

Components

Package	Items	Quantity
25 rxn	Annexin V-FAM	1 vial (0.125 mL)
50 rxn	Annexin V-FAM	1 vial (0.25 mL)
100 rxn	Annexin V-FAM	1 vial (0.5 mL)

Specifications

Reactivity

Mammalian

Application

Cell apoptosis assay

Detection Method

Fluorescent

Sample Type

Adherent cells | Suspension cells

Instruction

Tainan Headquarter

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 12 months under sterile conditions from date of receipt.

Protect from light.

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarter

☎ +886-6-2536677

✉ bd@leadgene.com.tw

Innovation & Research Center

☎ +886-2-27065528

CLD Center

☎ +886-6-2536677