

## NTP Set, 100 mM Solutions

Catalog Number	LDG0006RI
Package	1 mL*4 / Customized package

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

The NTP Set consists of 4 separate 100 mM solutions of ATP, UTP, CTP, and GTP (pH 7.5) as sodium salts.

#### Product Note

Mix thoroughly before use and dilute with nuclease-free water if necessary.

### Specifications

#### Application

In vitro transcription, RNA amplification, miRNA and siRNA synthesis

#### Concentration

100 mM

#### Purity

≥98% purity of each component confirmed by HPLC.

#### Form

Liquid. Clear aqueous solutions (pH 7.4).

### Instruction

#### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

#### Stability & Storage

This product is stable after storage at:

- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

#### Tainan Headquarter

+886-6-2536677

✉ [bd@leadgene.com.tw](mailto:bd@leadgene.com.tw)




#### Innovation & Research Center

+886-2-27065528

#### CLD Center

+886-6-2536677

**Disclaimer :** For Research Use or Further Manufacturing Only.

**Tainan Headquarter** +886-6-2536677 [bd@leadgene.com.tw](mailto:bd@leadgene.com.tw)**Innovation & Research Center** +886-2-27065528**CLD Center** +886-6-2536677