

2X Colorimetric RT-LAMP Master Mix

Catalog Number	LDG0020RF
Package	1 mL / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

Leadgene® 2X Colorimetric RT-LAMP Master Mix is an optimized master mix for reverse-transcription loop-mediated isothermal amplification (RT-LAMP) reactions. This product is a dual enzyme system, providing a rapid and sensitive detection in one pot. The amplified products can be visually observed by a metal indicator for calcium, hydroxynaphthol blue (HNB). Positive reactions will appear a color change from violet to sky blue, while negative reactions won't.

Product Note

- Primer concentration

Primer concentration can be titrated between 0.25X - 1X if undesired background signal appeared.

Specifications

Application

Loop-mediated Isothermal Amplification (LAMP)

Form

Liquid

Instruction

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

Stability & Storage

This product is stable after storage at:

- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Tainan Headquarter

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarter

☎ +886-6-2536677

✉ bd@leadgene.com.tw

Innovation & Research Center

☎ +886-2-27065528

CLD Center

☎ +886-6-2536677