

EffiStart™ 2X Probe qPCR Master Mix

Catalog Number L

Package

LDG0029RF

1 mL / Customized package

For full product information, images and publications, please visit our website.



Overview

Description

EffiStart[™] 2X Probe qPCR Master Mix is an optimized reaction mix for real-time polymerase chain reaction (qPCR) kit. This product contains Hot Start Taq DNA Polymerase and is suitable for probe-based detection and formulated as a 2-fold premix. Reaction can be simply set up by adding the DNA template, primers, and probes. This master mix does not contain ROX reference dye; it offers great convenience and minimizes the risk of cross-contamination.

Product Note

• Primer/Probe concentration

Final concentrations of 400 nM (each primer) are suitable for most reactions. To obtain optimal condition, primer concentration can be titrated between 0.2-1 µM.

A final concentration of 200 nM (probe) is suitable for most reactions. To obtain optimal condition, probe concentration can be titrated between 0.1-0.3 μ M.

• Annealing/Extension optimization

To obtain optimal condition, annealing/extension temperature can be adjusted between 55°C-65°C, annealing/extension time can be extended up to 60 sec.

Target length

Appropriate amplicon length should be arranged between 80-200 bp.

Specifications

Application

☑ bd@leadgene.com.tw

Reverse transcription, RT-PCR, qPCR

Form

Liquid

Tainan Headquarter Innovation & Research Center CLD Center © +886-6-2536677 © +886-2-27065528 © +886-6-2536677



Instruction

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

Stability & Storage

This product is stable after storage at:

 -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarter

Innovation & Research Center

CLD Center

& +886-6-2536677

& +886-2-27065528

& +886-6-2536677

☑ bd@leadgene.com.tw