

## HyFect® Transfection Reagent

**Catalog Number** LDG0001RA

**Package** 500 µL / Customized package

 Publications ( 27 )

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

HyFect® DNA Transfection Reagent is a low cell toxicity and non-liposome reagent optimized for DNA delivery. It provides effectively, reproducibly, and affordable benefits for scientific research.

### Specifications

#### Application

Cell transfection

#### Sample Type

Adherent cells | Suspension cells

### Instruction

#### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

#### Stability & Storage


This product is stable after storage at:

- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.


**Disclaimer :** For Research Use or Further Manufacturing Only.

#### Tainan Headquarter

 +886-6-2536677  +886-6-2531536

 [info@leadgene.com.tw](mailto:info@leadgene.com.tw)

#### Taipei Office

 +886-2-27065528