

HyLink™ HRP Conjugation Kit, 10 µg*2

Catalog Number	LDG0001RC
Package	10 µg*2 / Customized package

“ Publications (1)

For full product information, images and publications, please visit [our website](#).



Overview

Description

Horseradish peroxidase (HRP) is a 44 kDa glycoprotein that can be modified to crosslink with antibodies or other proteins, without compromising the enzymatic activity. HRP conjugates can be widely used in ELISA, western blotting and IHC due to its enzymatic characteristics. Leadgene HyLink™ HRP Conjugation Kit is designed for HRP conjugation of a small quantity (10 µg-1 mg) of antibody or protein. It provides a rapid and easy process with high efficiency to conjugate HRP to antibody or protein. The kit provides all the necessary components and high-activity HRP for conjugation.

Components

Package	Items	Quantity
10 µg × 2	HRP Conjugate	2 vials
	10× Modifier	1 vial
	Reduction Reagent	1 vial
	10× Quencher	1 vial

Instruction

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

Stability & Storage

This product is stable after storage at:

- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Tainan Headquarters

+886-6-2536677

✉ bd@leadgene.com.tw



Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarters +886-6-2536677 bd@leadgene.com.tw**Innovation & Research Center** +886-2-27065528**CLD Center** +886-6-2536677