

## Normal Mouse IgG

<b>Catalog Number</b>	LDG0001YD
<b>Package</b>	5 mg / Customized package

“ Publications ( 3 )

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

Mouse IgG is purified and isolated from pooled normal mouse serum by protein A. The purified normal mouse IgG may be used as an isotype control in various immunoassays, such as ELISA, WB, IP, FACS, ID, and IE.

### Specifications

#### Host

Mouse

#### Clonality

Polyclonal

#### Application

ELISA, WB, IP, FACS, ID, IE

#### Conjugation

Unconjugated

#### Concentration

Lot-dependent, please refer to product label.

#### Storage Buffer

Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.

#### Form

Liquid

### Instruction

#### Tainan Headquarters

+886-6-2536677

bd@leadgene.com.tw

#### Innovation & Research Center

+886-2-27065528

#### CLD Center

+886-6-2536677

---

### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

---

### Stability & Storage

This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

---

---

**Disclaimer :** For Research Use or Further Manufacturing Only.

---

#### Tainan Headquarters

☎ +886-6-2536677

✉ [bd@leadgene.com.tw](mailto:bd@leadgene.com.tw)

---

#### Innovation & Research Center

☎ +886-2-27065528

---

#### CLD Center

☎ +886-6-2536677