

Goat Anti-Rabbit IgG (H+L) Antibody

Catalog Number	LDG0003YE
Package	1 mg / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

Goat Anti-Rabbit IgG (H+L) Antibody is purified from goat, specifically react with rabbit IgG, heavy and light chain.

Specifications

Host

Goat

Clonality

Polyclonal

Immunogen

Purified rabbit IgG

Reactivity

Rabbit (Minor cross-reactivity to human in higher concentration)

Conjugation

Unconjugated

Concentration

Lot-dependent, please refer to product label.

Storage Buffer

Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.

Specificity

IgG

Form

Liquid

Instruction

Tainan Headquarters

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarters

☎ +886-6-2536677

✉ bd@leadgene.com.tw

Innovation & Research Center

☎ +886-2-27065528

CLD Center

☎ +886-6-2536677