

## Mouse IgG2a Kappa Isotype Control Antibody

Catalog Number	LDG0004YC
Package	Customized package / 1mL

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

Mouse IgG2a Kappa isotype control antibody is used to estimate non-specific binding of target primary antibodies.

### Specifications

#### Host

Mouse

#### Clonality

Monoclonal

#### Isotype

IgG2a

#### Application

ELISA, WB, IP, FACS, ID, IE

#### Conjugation

Unconjugated

#### Concentration

1 mg/mL

#### Storage Buffer

Phosphate Buffered Saline containing 0.03% ProClin 300, pH 7.4.

#### Form

Liquid

### Instruction

#### Tainan Headquarters

+886-6-2536677

bd@leadgene.com.tw

#### Innovation & Research Center

+886-2-27065528

#### CLD Center

+886-6-2536677

### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

### Stability & Storage

This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

**Disclaimer :** For Research Use or Further Manufacturing Only.

#### Tainan Headquarters

+886-6-2536677

bd@leadgene.com.tw

#### Innovation & Research Center

+886-2-27065528

#### CLD Center

+886-6-2536677