

## Chicken Anti-Sheep IgG (H+L) Antibody

Catalog Number	LDG0006YE
Package	1 mg / Customized package

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

Chicken Anti-Sheep IgG (H+L) Antibody is purified from chicken, specifically react with sheep IgG, heavy and light chain.

### Specifications

#### Host

Chicken

#### Clonality

Polyclonal

#### Immunogen

Purified sheep IgG

#### Reactivity

Sheep

#### Conjugation

Unconjugated

#### Concentration

Lot-dependent, please refer to product label.

#### Storage Buffer

Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.

#### Specificity

IgG

#### Form

Liquid

### Instruction

#### Tainan Headquarters

+886-6-2536677

bd@leadgene.com.tw

#### Innovation & Research Center

+886-2-27065528

#### CLD Center

+886-6-2536677

---

### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

---

### Stability & Storage

This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

---

---

**Disclaimer :** For Research Use or Further Manufacturing Only.

---

#### Tainan Headquarters

☎ +886-6-2536677

✉ [bd@leadgene.com.tw](mailto:bd@leadgene.com.tw)

---

#### Innovation & Research Center

☎ +886-2-27065528

---

#### CLD Center

☎ +886-6-2536677