

## Goat Anti-Rabbit IgG (H+L) Antibody (HRP)

**Catalog Number** LDG0014YE

**Package** 1 mL / Customized package

 Publications ( 32 )

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

Leadgene Goat Anti-Rabbit IgG (H+L) Antibody (HRP) is conjugated with high-activity horseradish peroxidase (HRP).

#### Product Note

Recommended dilution factor:

ELISA: 1:5000-25000

IHC: 1:500-2500

WB: 1:5000-25000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

### Specifications

#### Host

Goat

#### Clonality

Polyclonal

#### Immunogen

Purified rabbit IgG

#### Reactivity

Rabbit (Minor cross-reactivity to human in higher concentration)

#### Application

ELISA, IHC, WB

#### Conjugation


HRP

#### Tainan Headquarter

 +886-6-2536677  +886-6-2531536

 [info@leadgene.com.tw](mailto:info@leadgene.com.tw)

#### Taipei Office

 +886-2-27065528

**Storage Buffer**

Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300 and 50% glycerol, pH 7.4.

**Form**

Liquid

**Specificity**

IgG

**Instruction****Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

**Stability & Storage**

This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles and protect from light.

Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

**Disclaimer :** For Research Use or Further Manufacturing Only.

**Tainan Headquarter**

+886-6-2536677 +886-6-2531536

info@leadgene.com.tw

**Taipei Office**

+886-2-27065528