

Goat Anti-Rabbit IgG (H+L) Antibody (HRP)

Catalog Number	LDG0014YE
Package	1 mL / Customized package

“ Publications (32)

For full product information, images and publications, please visit [our website](#).



Overview

Description

Leadgene Goat Anti-Rabbit IgG (H+L) Antibody (HRP) is conjugated with high-activity horseradish peroxidase (HRP).

Product Note

Recommended dilution factor:

ELISA: 1:5000-25000

IHC: 1:500-2500

WB: 1:5000-25000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

Specifications

Host

Goat

Clonality

Polyclonal

Immunogen

Purified rabbit IgG

Reactivity

Rabbit (Minor cross-reactivity to human in higher concentration)

Application

ELISA, IHC, WB

Conjugation

HRP

Tainan Headquarters

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677

Storage Buffer

Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300 and 50% glycerol, pH 7.4.

Form

Liquid

Specificity

IgG

Instruction**Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles and protect from light.

Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarters

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677