

Goat Anti-Rabbit IgG Antibody Sepharose

Catalog Number LDG0024RD

Package 0.25 mL / 2 mL / 5 mL / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

Affinity purified antibodies were conjugated to NHS-sepharose. Leadgene Goat anti-rabbit IgG sepharose has 6-10 mg antibodies conjugated per mL of settled sepharose.

Product Note

The optimal protein concentration of the cell lysate should be determined by the end user.

Specifications

Host

Goat

Clonality

Polyclonal

Immunogen

Purified rabbit IgG

Reactivity

rabbit

Application

Immunoprecipitation

Conjugation

Sepharose

Specificity

Recognizes rabbit IgG

Form

Gel Immobilized 50% slurry in phosphate-buffered saline containing 0.05% sodium azide, pH 7.4.

Tainan Headquarter

+886-6-2536677 +886-6-2531536

info@leadgene.com.tw

Taipei Office

+886-2-27065528

Binding capacity

The binding and elution capacity of 1 mL sepharose are commonly more than 1 mg of IgG protein.

Instruction

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 12 months under sterile conditions from date of receipt.

Do not freeze.

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarter

☎ +886-6-2536677 📠 +886-6-2531536

✉ info@leadgene.com.tw

Taipei Office

☎ +886-2-27065528