

Sheep Anti-Chicken IgY (H+L) Antibody (TAMRA)

Catalog Number	LDG0049YE
Package	500 µL / 1 mL / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

TAMARA is a fluorescent dye for DNA, peptides and proteins. Leadgene Sheep Anti-Chicken IgY (H+L) Antibody (TAMRA) is conjugated with high-activity TAMARA

Product Note

Recommended dilution factor:

IFA: 1:200-1000

FACS: 1:200-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

Specifications

Host

Sheep

Clonality

Polyclonal

Immunogen

Purified chicken IgY

Reactivity

Chicken

Application

IFA, FACS

Conjugation

5(6)-TAMRA

Concentration

0.5 mg/mL

Storage Buffer

Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300, pH 7.4.

Tainan Headquarters

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677

Specificity

IgY

Form

Liquid

Instruction**Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 6 months under sterile conditions from date of receipt.

Do not freeze and protect from light.

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarters +886-6-2536677 bd@leadgene.com.tw**Innovation & Research Center** +886-2-27065528**CLD Center** +886-6-2536677