

## Goat Anti-Chicken IgY Antibody

Catalog Number	LDG0052YE
Package	1 mg / Customized package

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

Goat anti-chicken IgY antibody is purified from goats and specifically reacts with chicken IgY antibodies.

### Specifications

#### Host

Goat

#### Clonality

Polyclonal

#### Immunogen

Purified chicken IgY

#### Conjugation

Unconjugated

#### Concentration

Lot-dependent, please refer to product label.

#### Storage Buffer

Phosphate Buffered Saline containing 0.03% ProClin300, pH 7.4.

#### Specificity

IgY

#### Form

Liquid

### Instruction

#### Tainan Headquarters

+886-6-2536677

bd@leadgene.com.tw

#### Innovation & Research Center

+886-2-27065528

#### CLD Center

+886-6-2536677

### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

### Stability & Storage


This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.


Avoid repeated freeze/thaw cycles.

Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

## Image



ELISA titration of goat anti-chicken IgY antibody



SDS-PAGE of purified antibody under reduced condition.

Lane M: Protein marker

Lane 1: Goat anti-Chicken IgY (H+L) Antibody (2 µg, reduced)

**Disclaimer :** For Research Use or Further Manufacturing Only.