

## Anti-PGA5 Antibody [Clone 85-5]

**Catalog Number** LDG0054YA

**Package** 100 µg / Customized package

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

Pepsinogens are members of the peptidases that are secreted by gastric chief cells and are inactive polypeptide proenzymes known as zymogens. They are converted by gastric acid to active pepsins, which are involved in the digestion of proteins. Pepsinogen A5 (PGA5), an isozymogen of pepsinogen A, is the predominant endopeptidase in the gastric acid and is associated with atrophic gastritis and gastritis.

#### Product Note

Recommended dilution factor:

ELISA: 1:5000-20000

WB: 1:1000-5000

IFA: 1:100-1000

FACS: Assay dependent

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

### Specifications

#### Host

Mouse

#### Clonality

Monoclonal

#### Isotype

IgG1

#### Clone Name

clone 85-5

#### Tainan Headquarter

+886-6-2536677 +886-6-2531536

info@leadgene.com.tw

#### Taipei Office

+886-2-27065528

**Immunogen**

The immunogen information used to generate this antibody is proprietary.

**Application**

ELISA, WB, IFA, FACS

**Concentration**

1 mg/mL

**Specificity**

PGA5

**Reactivity**

Human

**Conjugation**

Unconjugated

**Storage Buffer**

Phosphate Buffered Saline containing 0.03% ProClin 300, pH 7.4.

**Form**

Liquid

**Instruction****Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

**Stability & Storage**

This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Suggestion: Divide antibody into several vials.

Keep only vials for usage at 2-8°C.


**Image****Tainan Headquarter**

+886-6-2536677 +886-6-2531536


info@leadgene.com.tw

**Taipei Office**


+886-2-27065528




Immunofluorescence analysis of Anti-PGA5 Antibody [Clone 85-5] HeLa cells were fixed in 4% PFA, permeabilized with PBS containing 0.1% Triton X-100. Cells were stained with mouse anti-PGA5 monoclonal antibody (1:100) followed by secondary antibodies (goat anti-Mouse IgG-iFluor 488, 1:200, green) and cell nuclei were stained with Hoechst 33342 (Blue).



FACS analysis of Anti-PGA5 Antibody [Clone 85-5] K562 cells were stained with mouse anti-PGA5 monoclonal antibody at 2 µg/ml (red) and without antibody control (black).



ELISA analysis of Anti-PGA5 Antibody [Clone 85-5] The specificity of anti-PGA5 antibody recognize recombinant proteins in ELISA. Red: Porcine pepsin; Orange: PGA5; Green: Pepsinogen I; Blue: Pepsinogen II; Black: BSA (negative control).



Western blotting analysis of Anti-PGA5 Antibody [Clone 85-5]  
 Recombinant proteins and mouse stomach lysates were stained with mouse anti-PGA5 monoclonal antibody at 1:5000 dilution. Lane 1: Purified human pepsinogen I (25 ng). Lane 2: Purified human pepsinogen II (25 ng). Lane 3: Mouse stomach lysate (30 µg). Lane 4: Human gastric juice (10 µL).

**Disclaimer :** For Research Use or Further Manufacturing Only.