

## Anti-SARS CoV-2 NP-Gold Conjugates

**Catalog Number** LDG0097YA

**Package** bulk

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

Mouse anti-SARS-CoV & CoV-2 NP mAb is conjugated with colloidal gold.

#### Product Note

Colloidal gold conjugate with Mouse anti-SARS-CoV & CoV-2 NP mAb | clone 102-7 (cat. LDG0083YA) and clone 88 (cat. LDG0096YA)

### Specifications

#### Reactivity

SARS-CoV & CoV-2

#### Application

LFIA

#### Conjugation

Colloidal gold

#### Concentration

OD520 ~ 0.75-0.80 (after 5X dilution)

#### Storage Buffer

Liquid / cherry red. Nano Gold-40 Conjugate Diluent, pH 9.0.

#### Specificity

Nucleocapsid protein

#### Form

Liquid

### Instruction

#### Tainan Headquarters

+886-6-2536677

bd@leadgene.com.tw

#### Innovation & Research Center

+886-2-27065528

#### CLD Center

+886-6-2536677

---

### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C.

---

---

### Stability & Storage

This product is stable after storage at:

- 2-8°C for 12 months under sterile conditions from date of receipt.

Briefly centrifuge vials before opening.

---

**Disclaimer :** For Research Use or Further Manufacturing Only.

---

#### Tainan Headquarters

☎ +886-6-2536677

✉ [bd@leadgene.com.tw](mailto:bd@leadgene.com.tw)

---

#### Innovation & Research Center

☎ +886-2-27065528

---

#### CLD Center

☎ +886-6-2536677