

Anti-Influenza A Virus NP Antibody [Clone IA04]

Catalog Number	LDG0153YA
Package	100 µg / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

Anti-Influenza A virus NP Antibody only recognizes nucleocapsid protein (NP) of influenza A viruses but not influenza B viruses. Influenza A viruses are RNA viruses, and their subtypes are labeled according to an H number (for the type of hemagglutinin) and an N number (for the type of neuraminidase).

Product Note

Recognize Influenza A virus NP in Lateral Flow, when clone IA04 antibody was paired with Anti-Influenza A virus NP Antibody [clone IA01] (cat. LDG0143YA).

Recommended dilution factor:
ELISA: 1:5000-20000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

Specifications

Host

Mouse

Clonality

Monoclonal

Isotype

IgG2a

Clone Name

clone IA04

Immunogen

Recombinant Influenza A virus nucleocapsid protein

Reactivity

Influenza A virus

Tainan Headquarters

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677

Application

ELISA, IFA, CLIA, LFIA

Concentration

1 mg/mL

Specificity

Nucleocapsid protein

Conjugation

Unconjugated

Storage Buffer

Phosphate Buffered Saline containing 0.2% ProClin 300, pH 7.4.

Form

Liquid

Instruction**Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.



Stability & Storage


This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.


Avoid repeated freeze/thaw cycles.

Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

Image**Tainan Headquarters** +886-6-2536677 bd@leadgene.com.tw**Innovation & Research Center** +886-2-27065528**CLD Center** +886-6-2536677



Immunofluorescence analysis of Anti-Influenza A virus NP Antibody [clone IA04] (1:500)



ELISA assay of Anti-Influenza A virus NP Antibody [clone IA04]

Disclaimer : For Research Use or Further Manufacturing Only.