

Anti-African Swine Fever Virus (ASFV) p30 Antibody [Clone 16E1]

Catalog Number LDG0161YA

Package 100 µg / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

Anti-African Swine Fever Virus (ASFV) p30 Antibody [Clone 16E11] recognizes ASFV p30 proteins. African swine fever (ASF) is a contagious septicemic viral disease of pigs. The pathogen of this disease is ASFV, a large DNA virus from the Asfarviridae family. The ASFV p30 protein is part of the inner membrane of the virion. Therefore, p30 protein is an important target for serology and antigen detection in ASF.

Product Note

Recognize African Swine Fever Virus p30 in ELISA, when monoclonal antibody was paired with Anti-African Swine Fever Virus (ASFV) p30 Antibody [Clone 18H1] (cat. LDG0162YA).

Recommended dilution factor:

ELISA: 1:5000-20000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

Specifications

Host

Mouse

Clonality

Monoclonal

Isotype

IgG1

Clone Name

clone 16E1

Immunogen

Recombinant African Swine Fever Virus p30 protein

Reactivity

African swine fever virus

Tainan Headquarters

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677

Application

ELISA, CLIA

Concentration

Lot-dependent, please refer to product label.

Specificity

p30 protein

Conjugation

Unconjugated

Storage Buffer

Phosphate Buffered Saline containing 0.03% ProClin 300, pH 7.4.

Form

Liquid

Instruction**Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

Stability & Storage


This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarters +886-6-2536677 bd@leadgene.com.tw**Innovation & Research Center** +886-2-27065528**CLD Center** +886-6-2536677