

## Mouse IL-36RA, His Tag, E. coli

Catalog Number	LDG059PME
Package	5 µg / 20 µg / 100 µg / Customized package

For full product information, images and publications, please visit [our website](#).



### Specifications

#### Species of Origin

Mouse

#### Affinity Tag

His Tag (C-term)

#### Purity

>98% as determined by SDS-PAGE analysis.

#### Activity

Measure by its ability to inhibit IL-36 gamma-induced IL-6 secretion in 3T3 cells. The ED<sub>50</sub> for this effect is <2 µg/mL.

#### Form

Lyophilized

#### Expression System

Escherichia coli

#### Storage Buffer

Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.

#### Molecular weight

The protein has a calculated MW of 17.81 kDa. The protein migrates as 17-25 kDa under reducing condition (SDS-PAGE analysis).

#### Endotoxin Level

<0.1 EU per 1 µg of the protein by the LAL method.

### Background

#### Tainan Headquarters

+886-6-2536677

bd@leadgene.com.tw

#### Innovation & Research Center

+886-2-27065528

#### CLD Center

+886-6-2536677

### Background

Interleukin 36 receptor antagonist (IL-36RA) is a 17.21 kDa cytokine with 155 amino acid residues and a member of the IL-36 family. It is mainly secreted from the stomach, lungs, and kidneys. IL-36RA binds to the interleukin-36 receptor and plays a critical role in the anti-inflammatory effect by mediating MAPK and NF- $\kappa$ B signaling pathways. IL-36RA inhibits proinflammatory effect via suppressing production of IL-6, IL-17, and interferon-gamma and promoting production of IL-4.

### Uniprot ID

#Q9QYY1

### Synonyms

Interleukin-36 receptor antagonist protein, IL-36Ra, Interleukin-1 HY1, IL-1HY1, Interleukin-1 delta, IL-1 delta, Interleukin-1 family member 5, IL-1F5, Interleukin-1 homolog 3, IL-1H3, Interleukin-1-like protein 1, IL-1L1

### Sequence Note

Met2-Asp156

## Instruction

### Reconstitution

It is recommended to reconstitute the lyophilized protein in sterile H<sub>2</sub>O to a concentration not less than 200  $\mu$ g/mL and incubate the stock solution for at least 20 min to ensure sufficient re-dissolved.

### Stability & Storage

This product is stable after storage at:

- -20°C for 12 months in lyophilized state from date of receipt.
- -20°C or -80°C for 1 month under sterile conditions after reconstitution.

Avoid repeated freeze/thaw cycles.

### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

## Image

### Tainan Headquarters

+886-6-2536677


bd@leadgene.com.tw

### Innovation & Research Center

+886-2-27065528

### CLD Center

+886-6-2536677



SDS-PAGE analysis of recombinant mouse IL-36RA.

**Disclaimer :** For Research Use or Further Manufacturing Only.