

Human ENO1, His Tag, E. coli

Catalog Number	LDG174PHE
Package	5 µg / 20 µg / 100 µg / Customized package

“ Publications (2)

For full product information, images and publications, please visit [our website](#).



Specifications

Species of Origin

Human

Expression System

Escherichia coli

Affinity Tag

His Tag (N-term)

Storage Buffer

Lyophilized from a 0.2 µm filtered solution of PBS with 7 mM MgSO₄, pH 7.2

Purity

>98% as determined by SDS-PAGE analysis.

Molecular weight

The protein has a calculated MW of 47.85 kDa.
The protein migrates as 55 kDa under reducing condition (SDS-PAGE analysis).

Endotoxin Level

<0.1 EU per 1 µg of the protein by the LAL method.

Form

Lyophilized

Background

Tainan Headquarters

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677

Background

Enolase 1 (ENO1) also named alpha-enolase which is a glycolytic enzyme. It could convert the NAD⁺ to NADH by producing two ATP. ENO1 is a homodimer which is composed by two isozymes of enolase (2 α , 2 ν , or 2 β). It also plays an important role in cancer model which can promote tumor cell proliferation and migration by the PI3K signaling pathway. ENO1 is a biomarker of prognostic and diagnostic cancer.

Uniprot ID

P06733 2

Synonyms

Alpha-enolase, 2-phospho-D-glycerate hydro-lyase, C-myc promoter-binding protein, Phosphopyruvate hydratase, Plasminogen-binding protein

Sequence Note

Ser2-Lys434

Instruction

Reconstitution

It is recommended to reconstitute the lyophilized protein in sterile H₂O to a concentration of 200 μ g/mL and incubate the stock solution for at least 20 min to ensure sufficient re-dissolved.

Stability & Storage

This product is stable after storage at:

- -20°C for 12 months in lyophilized state from date of receipt.
- -20°C or -80°C for 1 month under sterile conditions after reconstitution.

Avoid repeated freeze/thaw cycles.

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

Image

Tainan Headquarters

+886-6-2536677


bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677



SDS-PAGE analysis of recombinant human ENO1.

Disclaimer : For Research Use or Further Manufacturing Only.