

Isotype Control List

Product Name	Buffer	Endotoxin	Purity	Package
Anti-HEL Rat IgG1-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Rabbit IgG-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Mouse IgG2C-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Mouse IgG2B-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL scFv mouse IgG2A-Fc(L234A/L235A) Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Mouse IgG2A(L234A/L235A)-Kappa Isotype-control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Mouse IgG2A-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Mouse IgG1-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Mouse IgG1(D265A)-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG4(S228P)-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG4(S228P/L235E)-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG4-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG3-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG2-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG1-Lambda1 Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG1(S239D/I332E)-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG1(L234A /L235A)-Kappa Isotype-control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG1(L234A /L235A/P329G)-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG1(N297A)-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG1(K214R)-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG1 CH1-Human IgG2 Fc-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL HlgG4(S228P/R409K)-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG1-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG1(L234F/L235E/P331S)-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG1(E356D/M358L)-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG1(K214R/K447A)-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG1(Fab-His) Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
anti-HEL Human IgG4(S228P/F234A/L235A)-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG4(S228P/F234A/L235A)-Lambda1 Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL Human IgG1(S239D/A330L/I332E))-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(N297S)-Kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(N297Q)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(P331S/E430G)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(E233A/L235A)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(K214R/N297A)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(D265A/N297A)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(K214R/L234A/L235A)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(K409R/E356D/M358L)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(N297A/E356D/M358L)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(P329G/E356D/M358L)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(L234F/L235E)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(D265A)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(K214R/L234F/L235E/P331S)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(K214R/E356D/M358L)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(L234A/L235A/G237A)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(L234A/L235A/F405L)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(D270N)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(F405L)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(E356D/M358L/L234A/L235A)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG1(E356D/M358L/L234A/L235A/F405L)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG4(S228P/L235E/P329G)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG4(S228P/L235E/D265A/R409K)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg

Product Name	Buffer	Endotoxin	Purity	Package
Anti-HEL-Human IgG4(S228P/F234A/L235A/R409K)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG4(S228P/L235E/R409K)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Human IgG4(S228P/D265A/L309V/R409K)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Mouse IgG1(N297A)-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL-Mouse IgG3-kappa Isotype control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL VHH-human IgG1 Fc Isotype Control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL VHH Isotype Control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL VHH-His Isotype Control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL VHH- Mouse IgG1 Fc Isotype Control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg
Anti-HEL VHH- Mouse IgG2A Fc Isotype Control	PBS, pH 7.4	<1EU/mg	>95%	1mg/5mg/25mg/50mg/100mg