

## Normal Rabbit IgG

**Catalog Number** LDG0002YD

**Package** 50 µg / 5 mg / 10 mg / Customized package

 Publications ( 3 )

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

Rabbit IgG is isolated and purified from pooled normal rabbit serum by protein A. The purified normal rabbit IgG may be used as an isotype control in various immunoassays, such as ELISA, WB, IP, FACS, ID, and IE.

### Specifications

#### Host

Rabbit

#### Clonality

Polyclonal

#### Application

ELISA, WB, IP, FACS, ID, IE

#### Conjugation

Unconjugated

#### Concentration

Lot-dependent, please refer to product label.

#### Buffer


Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.

#### Form

Liquid


### Instruction

#### Tainan Headquarter

 +886-6-2536677  +886-6-2531536

 [info@leadgene.com.tw](mailto:info@leadgene.com.tw)

#### Taipei Office

 +886-2-27065528

---

### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

---

### Stability & Storage

"This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Suggestion: Divide antibody into several vials.

Keep only vials for usage at 2-8°C."

---

**Disclaimer :** For Research Use or Further Manufacturing Only.

---

### Tainan Headquarter

☎ +886-6-2536677 📠 +886-6-2531536

✉ [info@leadgene.com.tw](mailto:info@leadgene.com.tw)

---

### Taipei Office

☎ +886-2-27065528