

## Goat Anti-Rabbit IgG (H+L) Antibody

<b>Catalog Number</b>	LDG0003YE
<b>Package</b>	1 mg / Customized package

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

Goat Anti-Rabbit IgG (H+L) Antibody is purified from goat, specifically react with rabbit IgG, heavy and light chain.

### Specifications

#### Host

Goat

#### Clonality

Polyclonal

#### Immunogen

Purified rabbit IgG

#### Reactivity

Rabbit (Minor cross-reactivity to human in higher concentration)

#### Conjugation

Unconjugated

#### Concentration

Lot-dependent, please refer to product label.

#### Buffer

Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.

#### Specificity

IgG

#### Form

Liquid

### Instruction

#### Tainan Headquarter

+886-6-2536677

bd@leadgene.com.tw

#### Innovation & Research Center

+886-2-27065528

#### CLD Center

+886-6-2536677

---

### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

---

### Stability & Storage

This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

---

---

**Disclaimer :** For Research Use or Further Manufacturing Only.