

CTP Solution (100 mM)

Catalog Number	LDG0004RI
Package	1 mL / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

Leadgene Cytidine-5'-triphosphate (CTP) is supplied as a 100 mM aqueous solution at pH 7.5.

Product Note

Mix thoroughly before use and dilute with nuclease-free water if necessary.

Specifications

Application

In vitro transcription, RNA amplification, miRNA and siRNA synthesis

Concentration

100 mM

Purity

≥98% purity of each component confirmed by HPLC.

Form

Liquid. Clear aqueous solutions (pH 7.4).

Instruction

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

Stability & Storage

This product is stable after storage at:

- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Tainan Headquarter

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarter

☎ +886-6-2536677

✉ bd@leadgene.com.tw

Innovation & Research Center

☎ +886-2-27065528

CLD Center

☎ +886-6-2536677