

## Normal Chicken IgY

<b>Catalog Number</b>	LDG0004YD
<b>Package</b>	1 mg / 5 mg / Customized package

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

Chicken IgY is isolated and purified from pooled normal chicken serum by protein A. The purified normal chicken IgY may be used as an isotype control in various immunoassays, such as ELISA, WB, IP, FACS, ID, and IE.

### Specifications

#### Host

Chicken

#### Clonality

Polyclonal

#### Application

ELISA, WB, IP, FACS, ID, IE

#### Conjugation

Unconjugated

#### Concentration

Lot-dependent, please refer to product label.

#### Buffer

Phosphate Buffered Saline containing 0.03% Proclin 300, pH 7.4.

#### Form

Liquid

### Instruction

#### Tainan Headquarter

+886-6-2536677

bd@leadgene.com.tw

#### Innovation & Research Center

+886-2-27065528

#### CLD Center

+886-6-2536677

## Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

## Stability & Storage


"This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C."

## Image



SDS-PAGE of purified antibody under reduced condition.

Lane M: Protein marker

Lane 1: Normal Chicken IgY (2 µg, reduced)

**Disclaimer :** For Research Use or Further Manufacturing Only.