

Donkey Anti-Goat IgG (H+L) Antibody

Catalog Number	LDG0004YE
Package	1 mg / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

Donkey Anti-Goat IgG (H+L) Antibody is purified from donkey, specifically react with goat IgG, heavy and light chain.

Specifications

Host

Donkey

Clonality

Polyclonal

Immunogen

Purified goat IgG

Reactivity

Goat (Minor cross-reactivity to human in higher concentration)

Conjugation

Unconjugated

Concentration

1 mg/mL

Buffer

Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.

Specificity

IgG

Form

Liquid

Instruction

Tainan Headquarter

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

Disclaimer : For Research Use or Further Manufacturing Only.