

## Mouse IgG2b Kappa Isotype Control Antibody

**Catalog Number** LDG0005YC

**Package** 1mL / Customized package

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

Mouse IgG2b Kappa isotype control antibody is used to estimate non-specific binding of target primary antibodies.

### Specifications

#### Host

Mouse

#### Clonality

Monoclonal

#### Isotype

IgG2b

#### Application

ELISA, WB, IP, FACS, ID, IE

#### Conjugation

Unconjugated

#### Concentration

1 mg/mL

#### Buffer

Phosphate Buffered Saline containing 0.03% ProClin 300, pH 7.4.

#### Form

Liquid

### Instruction

#### Tainan Headquarter

+886-6-2536677 +886-6-2531536

info@leadgene.com.tw

#### Taipei Office

+886-2-27065528

---

### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

---

### Stability & Storage

This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Suggestion: Divide antibody into several vials.

Keep only vials for usage at 2-8°C.

---

**Disclaimer :** For Research Use or Further Manufacturing Only.

---

#### Tainan Headquarter

☎ +886-6-2536677    📠 +886-6-2531536

✉ [info@leadgene.com.tw](mailto:info@leadgene.com.tw)

---

#### Taipei Office

☎ +886-2-27065528