

Goat Anti-Human IgG (Fc Fragment) Antibody

Catalog Number LDG0011YE

Package 1 mg / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

Leadgene Goat anti-Human IgG, Fc Fragment is purified from goat antiserum by immunoaffinity chromatography. Minimal cross-reactivity with other isotype or species.

Specifications

Host

Goat

Clonality

Polyclonal

Immunogen

Gamma immunoglobulin Fc region

Reactivity

Human

Conjugation

Unconjugated

Concentration

1 mg/mL

Buffer

Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.

Specificity

IgG, Fc Fragment

Form

Liquid

Instruction

Tainan Headquarter

+886-6-2536677 +886-6-2531536

info@leadgene.com.tw

Taipei Office

+886-2-27065528

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Suggestion: Divide antibody into several vials.

Keep only vials for usage at 2-8°C.

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarter

☎ +886-6-2536677 📠 +886-6-2531536

✉ info@leadgene.com.tw

Taipei Office

☎ +886-2-27065528