

## LG-1 2X Diluent Buffer

**Catalog Number** LDG0017RH

**Package** 10 mL / 50 mL / 100 mL

“ Publications ( 1 )

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

LG-1 2X Diluent Buffer is used for dilution of biological samples or antibodies. LG-1 2X Diluent Buffer can reduce matrix effect and non-specific binding in various immunoassays, such as ELISA, WB and IHC.

### Specifications

#### Form

Liquid. LG-1 2X Diluent Buffer was aseptically filtered through a 0.22 µm filter.

### Instruction

#### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

#### Stability & Storage

This product is stable after storage at:

- 2-8°C for 12 months in lyophilized state from date of receipt.

**Disclaimer :** For Research Use or Further Manufacturing Only.

**Tainan Headquarter**

+886-6-2536677

bd@leadgene.com.tw

**Innovation & Research Center**

+886-2-27065528

**CLD Center**

+886-6-2536677