

## LG-2 2X Diluent Buffer

Catalog Number	LDG0018RH
Package	10 mL / 50 mL / 100 mL

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

LG-2 2X Diluent Buffer is used for dilution of biological samples or antibodies. LG-2 2X Diluent Buffer can reduce matrix effect and non-specific binding in various immunoassays, such as ELISA, WB and IHC.

### Specifications

#### Form

Liquid. LG-2 2X Diluent Buffer was aseptically filtered through a 0.22  $\mu$ m filter.

### Instruction

#### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

#### Stability & Storage

This product is stable after storage at:

- 2-8°C for 12 months in lyophilized state from date of receipt.

**Disclaimer :** For Research Use or Further Manufacturing Only.

#### Tainan Headquarter

+886-6-2536677

✉ [bd@leadgene.com.tw](mailto:bd@leadgene.com.tw)

#### Innovation & Research Center

+886-2-27065528

#### CLD Center

+886-6-2536677