

Sheep Anti-Chicken IgY (H+L) Antibody (HRP)

Catalog Number LDG0018YE

Package 1 mL / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

Leadgene Sheep Anti-Chicken IgY (H+L) Antibody (HRP) is conjugated with high-activity horseradish peroxidase (HRP).

Product Note

Recommended dilution factor:

ELISA: 1:5000-25000

IHC: 1:500-2500

WB: 1:5000-25000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

Specifications

Host

Sheep

Clonality

Polyclonal

Immunogen

Purified chicken IgY

Reactivity

Chicken

Application

ELISA, IHC, WB

Conjugation

HRP

Tainan Headquarter

+886-6-2536677 +886-6-2531536

info@leadgene.com.tw

Taipei Office

+886-2-27065528

Buffer

Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300 and 10% glycerol, pH 7.4.

Form

Liquid

Specificity

IgY

Instruction**Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles and protect from light.

Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarter

+886-6-2536677 +886-6-2531536

info@leadgene.com.tw

Taipei Office

+886-2-27065528