

LG-3 2X Diluent Buffer

Catalog Number LDG0019RH

Package 10 mL / 50 mL / 100 mL

For full product information, images and publications, please visit [our website](#).



Overview

Description

LG-3 2X Diluent Buffer is used for dilution of biological samples or antibodies. LG-3 2X Diluent Buffer can reduce matrix effect and non-specific binding in various immunoassays, such as ELISA, WB and IHC.

Specifications

Form

Liquid. LG-3 2X Diluent Buffer was aseptically filtered through a 0.22 µm filter.

Instruction

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 12 months in lyophilized state from date of receipt.

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarter

☎ +886-6-2536677 📠 +886-6-2531536

✉ info@leadgene.com.tw

Taipei Office

☎ +886-2-27065528