

Anti-mCherry Antibody Sepharose

Catalog Number	LDG0021RD
Package	0.25 mL / 2 mL / 5 mL / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

Anti-mCherry VHH antibody was conjugated to NHS-sepharose beads. Leadgene Mouse anti-mcherry mAb sepharose has 6-10 mg mAb conjugated per mL of settled sepharose.

Product Note

The optimal protein concentration of the cell lysate should be determined by the end user.

Specifications

Clonality

VHH

Isotype

IgG2a

Immunogen

Recombinant mCherry protein

Application

Immunoprecipitation, Immunoaffinity Purification

Conjugation

Sepharose

Specificity

Recognizes RFP, DsRed and mCherry tag on either the N- or C-terminus of the recombinant fusion protein.

Sample Type

Cell lysate

Form

Gel Immobilized 50% slurry in 20% ethanol.

Binding capacity

The binding and elution capacity of 1 mL sepharose are commonly more than 1 mg of RFP, mCherry or DsRed fusion protein.

Tainan Headquarter

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677

Instruction

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 12 months under sterile conditions from date of receipt.

Do not freeze.

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarter

+886-6-2536677

✉ bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677