

Goat Anti-Mouse IgG Antibody Sepharose

Catalog Number	LDG0023RD
Package	0.25 mL / 2 mL / 5 mL / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

Affinity purified antibodies were conjugated to NHS-sepharose beads. Leadgene Goat anti-Mouse IgG sepharose has 6-10 mg pAb conjugated per mL of settled sepharose.

Product Note

The optimal protein concentration of the cell lysate should be determined by the end user.

Specifications

Host

Goat

Clonality

Polyclonal

Immunogen

Purified mouse IgG

Reactivity

mouse

Application

Immunoprecipitation

Conjugation

Sepharose

Specificity

Recognizes mouse IgG

Form

Gel Immobilized 50% slurry in phosphate-buffered saline containing 0.05% sodium azide, pH 7.4.

Binding capacity

The binding and elution capacity of 1 mL sepharose are commonly more than 1 mg of IgG protein.

Tainan Headquarter

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677

Instruction

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 12 months under sterile conditions from date of receipt.

Do not freeze.

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarter

+886-6-2536677

✉ bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677