

Anti-SARS CoV-2 NP-Gold Conjugates

Catalog Number LDG0097YA

Package bulk

For full product information, images and publications, please visit [our website](#).



Overview

Description

Mouse anti-SARS-CoV & CoV-2 NP mAb is conjugated with colloidal gold.

Product Note

Colloidal gold conjugate with Mouse anti-SARS-CoV & CoV-2 NP mAb | clone 102-7 (cat. LDG0083YA) and clone 88 (cat. LDG0096YA)

Specifications

Reactivity

SARS-CoV & CoV-2

Application

LFIA

Conjugation

Colloidal gold

Concentration

OD520 ~ 0.75-0.80 (after 5X dilution)

Buffer

Liquid / cherry red. Nano Gold-40 Conjugate Diluent, pH 9.0.

Specificity

Nucleocapsid protein

Form

Liquid

Instruction

Tainan Headquarter

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 12 months under sterile conditions from date of receipt.

Briefly centrifuge vials before opening.

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarter

☎ +886-6-2536677

✉ bd@leadgene.com.tw

Innovation & Research Center

☎ +886-2-27065528

CLD Center

☎ +886-6-2536677